

Keeler Loupes

The brand choice for Dental, Surgical and Veterinary Professionals



Now
with NEW
Keeler
Frames




Keeler

Take a closer look

Magnification allows you to do a better job. If you can't see it, you can't work on it. The benefits of Keeler loupes are immediately clear: improved working posture, reduced eye strain and now, with the new K-L.E.D. head light, your working view is even clearer.

Laser aligned optics for magnification of superb optical quality, a lightweight aluminium assembly and a precise pupil alignment system make Keeler loupes the perfect choice for dental, surgical and veterinary professionals.



SuperVu Galilean system

Delivering 2.0x, 2.5x and 3.0x magnification.

The true benefit of this system is ease of use. Movement whilst maintaining focal clarity is possible due to superior depth of field. The lightweight aluminium bar with individual pupillary adjustment provides the mounted loupes with stability and optimal focal alignment, eliminating headache and eyestrain.



SuperVu Galilean system

SuperVu Hi-Res Galilean system NEW

Delivering 3.0x magnification.

Keeler have designed a new high resolution Galilean loupe which exceeds the standards previously attained. The 5 element loupe gives a true 3x magnification with a 15% increased field of view and a visually flat field from edge to edge.



NEW SuperVu
Hi-Res Galilean system



Protective Caps

SuperVu Galilean system Part Numbers

2.0x Magnification

2112-P-1002 SuperVu 2.0x @ 34cm working distance

2112-P-1003 SuperVu 2.0x @ 46cm working distance

2.5x Magnification

2112-P-1010 SuperVu 2.5x @ 34cm working distance

2112-P-1011 SuperVu 2.5x @ 42cm working distance

2112-P-1012 SuperVu 2.5x @ 46cm working distance

2112-P-1013 SuperVu 2.5x @ 50cm working distance

3.0x Magnification

2112-P-1022 SuperVu 3.0x @ 34cm working distance

2112-P-1023 SuperVu 3.0x @ 42cm working distance

2112-P-1024 SuperVu 3.0x @ 46cm working distance

2112-P-1025 SuperVu 3.0x @ 50cm working distance

SuperVu Hi-Res Galilean system Part Numbers

3.0x Magnification

2112-P-1050 SuperVu Hi Res 3.0x @ 34cm working distance

2112-P-1051 SuperVu Hi Res 3.0x @ 42cm working distance

2112-P-1052 SuperVu Hi Res 3.0x @ 46cm working distance

2112-P-1053 SuperVu Hi Res 3.0x @ 50cm working distance

SuperVu accessories

2112-P-8049 Protective Caps

2112-P-8057 +1D Near distance Caps

Let Keeler loupes take the strain

Keeler Loupes are comfortable to wear and work with, helping to reduce eyestrain. By bringing your patients closer, you can improve your posture and avoid stooping.

Patient comfort

Patients will feel more comfortable and at ease with increased personal space as you will work from a greater distance.

XL Advantage system

Delivering 3.5x, 4.5x and 5.5x magnification.

Higher magnification for advanced surgical procedures. Prismatic quality guarantees a true image. Edge to edge clarity and precise colour rendition is ideal for intricate applications. The XL's precision lenses are treated with an anti-reflective coating to produce a superior image.

Multiple magnification levels

Keeler Loupes offer performance matched to your requirements. Multiple magnification levels enable you to see more, spot problems earlier and offer even higher levels of precision treatment.



XL Advantage system

XL Advantage system Part Numbers

3.5x Magnification

2113-P-2004 XL Advantage 3.5x @ 34cm working distance

2113-P-2005 XL Advantage 3.5x @ 42cm working distance

2113-P-2006 XL Advantage 3.5x @ 46cm working distance

4.5x Magnification

2113-P-2008 XL Advantage 4.5x @ 34cm working distance

2113-P-2009 XL Advantage 4.5x @ 42cm working distance

2113-P-2010 XL Advantage 4.5x @ 46cm working distance

5.5x Magnification

2113-P-2012 XL Advantage 5.5x @ 34cm working distance

2113-P-2013 XL Advantage 5.5x @ 42cm working distance

2113-P-2014 XL Advantage 5.5x @ 46cm working distance

NEW Keeler Frames

All Keeler loupes are available on the new Keeler frame.

Stylish, comfortable and attractive, the new Keeler frame is available in a choice of bright colours as shown. The frames will fit over many styles of glasses, or with the free insert, can be easily adapted to all prescriptions.



Sky Blue



Black



Blue



Red



K-L.E.D. gives long lasting white light

The K-L.E.D. head light produces a homogeneous light patch of up to 18000 lux.

The 3W L.E.D. is mounted coaxially on the loupe bar and is powered by a lithium battery pack allowing 4 hours of continuous use. The K-L.E.D. is available in two packages; a twin charger with a tray to hold loupes when not in use that will charge two batteries simultaneously, or a kit bag with a compact single charger that gives the K-L.E.D. even greater portability. Both are available with either a universal clip for attachment to most brands of metal framed loupes or glasses, or a fitting for all current and previous model Keeler loupes.



NEW K-L.E.D. and SuperVu Galilean system



K-L.E.D. with single charger and bag



K-L.E.D. and SuperVu Galilean system in twin charger tray

Part Numbers

2004-P-1000	Keeler fit K-L.E.D. system with twin charger tray
2004-P-1003	Universal fit K-L.E.D. system with twin charger tray
2004-P-1005	Universal fit K-L.E.D. system with single charger and bag
2004-P-1006	Keeler fit K-L.E.D. system with single charger and bag



Designed for comfort

Choose from the new Keeler Sport frame or a headband mounted system to mount your loupes. Adaptable to all prescriptions.



Call FREEPHONE 0800 521251 to order or arrange your free demonstration today

Keeler Ltd.
Clewer Hill Road, Windsor
Berkshire, SL4 4AA
Tel: +44 (0) 1753 857177
Fax: +44 (0) 1753 827145
Email: info@keeler.co.uk
Website: www.keeler.co.uk

Keeler Scotland
25 Deerdykes View
Westfield Estate
Cumbernauld G68 9HN
Tel: +44 (0) 1236 721214
Fax: +44 (0) 1236 721231